

ABOUT RECENDE

Zhaoqing Gaoyao Recende Medical Equipment Co.,LTD is a professional manufacture of **Portable Dental Unit ,Dental Unit ,Dental Air Compressor, Dental X-ray Unit , Autoclave / Sterilizer** that is specialized in the development manufacture and sell of high quality Dental equipment. During these years , with the joint effort of all staff and supports of our customers

RECENDE has grown rapidly. We has advanced technology terms, user - friendly design, strict quality control & good after-sales service. **RECENDE** products have been sold to Asia, Africa, Oceania, South America, Europe and received good feedback from our agents and customers.

We are confident to build the most affordable high-tech products and keep you stay on the leading-edge of dental technology.

We will do our best for your health!



PORTABLE DENTAL UNIT



BL-604

BL-605

Army
Medical
Hospitals
Home visits
Veterinarian
Dental clinic
Visit of care home
School oral health
Humanitarian aid

Standard Equipment

Silent Oilless Air Compressor	1pc
Three Way Syringe	1pc
High Speed Handpiece Tube	1pc
Low Speed Handpiece Tube	1pc
Stainless Steel Air Tank 6L	1pc
Adjust Strong Suction System	1pc
Water Bottle	1pc
Drain Bottle	1pc
Foot Control	1pc
Built in Ultrasonic Scaler	1pc
Built in LED Light Curing Unit	1pc

Specification

Power Supply: AC110V or AC220V
Power: 550W
Pressure: 0.8Mpa
Noise: 38dB
G.W.: 33Kgs
Dimension: 53x40x73(cm)

Options

High Speed Handpiece
Low Speed Handpiece
Portable Dental Chair with light

PORTABLE DENTAL UNIT



BL-601 & BL-602



BL-603



Back size

Classification	Details	BL-601	BL-602	BL-603
Specification	Power Supply(AC)	110V/220V	110V/220V	110V/220V
	Power(W)	550	550	550
	Pressure(Mpa)	0.8	0.8	0.8
	Noise(dB)	39	40	40
	G.Weight(Kg)	19	21	22
	Dimension(CM)	48x40x46	48x40x46	47x37x42
Standard Equipment	High Speed Handpiece	★	★	★
	Low Speed Handpiece	★	★	★
	Three Way Syringe	★	★	★
	Silent Oilless Air Compressor	★	★	★
	Stainless Steel Air Tank	★	★	★
	Water Bottle	★	★	★
	Suction System	★	★	★
	Drain Bottle	★	★	★
	Foot Control	★	★	★
	Pull rod	★	★	★
Options	High Speed Handpiece	▲	▲	▲
	Low Speed Handpiece	▲	▲	▲
	Built in LED Light Curing Unit	▲	▲	▲
	Built in Ultrasonic Scaler	▲	▲	▲
	Portable Dental Chair with light	▲	▲	▲

PORTABLE DENTAL UNIT

Standard Equipment:

LED Inner Electrical Motor	1set
Built in LED Ultrasonic Scaler	1set
Three Way Syringe	1set
Silent Oilless Air Compressor	1set
Stainless Steel Water Tank	1set
Stainless Steel Air Tank	1set
Strong Suction System	1pc
Drain Bottle	1pc
Foot Control	1pc

Options:

Contra-angle Handpiece 1:1, 1:5, 16:1
Portable Dental Chair with light



Digital Display



BL-608

Portable Dental Chair



BL-612



BL-613



Bag for Portable Dental Chair

DENTAL UNIT

BL-801



Standard Equipment

Operation Cold Lamp
Foot control
Doctor stool
PU upholstery Dental Chair
24V noiseless DC motor
Water heater system

Tempered glass spittoon
Two three-way syringe(cold/hot water)
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

BL-801A



Standard Equipment

LED Operation Lamp
Foot control
Doctor stool
PU upholstery Dental Chair
24V noiseless DC motor
Water heater system

Fixation ceramic spittoon
Two three-way syringe(cold/hot water)
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

DENTAL UNIT

BL-802



Standard Equipment

Operation Cold Lamp
Foot control
Doctor stool
PU upholstery Dental Chair
24V noiseless DC motor
Water heater system

Tempered glass spittoon
Two three-way syringe(cold/hot water)
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

BL-802A



Standard Equipment

Operation Cold Lamp
Foot control
Doctor stool
PU upholstery Dental Chair
24V noiseless DC motor
Water heater system

High quality move spittoon
Two three-way syringe(cold/hot water)
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

DENTAL UNIT

BL-803



Standard Equipment

LED Operation Lamp
Foot control
Doctor stool
24V noiseless DC motor
Water heater system
Tempered glass spittoon

Imported air & water tube
Two three-way syringe(cold/hot water)
Seamless PU upholstery Dental Chair
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

BL-803A



Standard Equipment

LED Operation Lamp
Foot control
Doctor stool
24V noiseless DC motor
Water heater system
Tempered glass spittoon

Upper hanging tools
Two three-way syringe(cold/hot water)
Seamless PU upholstery Dental Chair
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

DENTAL UNIT

BL-804



Standard Equipment

LED Operation Lamp
Foot control
Doctor stool
Taiwan motor
Water heater system
Luxury ceramic spittoon
Low hanging tools

Imported air & water tube
Two three-way syringe(cold/hot water)
Seamless PU upholstery Dental Chair
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Two high speed and one low speed handpieces tubes

BL-805



Standard Equipment

Foot control
Doctor stool
Water heater system
Luxury glass spittoon
LED Operation Lamp with sensor
Imported air & water tube
Taiwan motor.

Two three-way syringe(cold/hot water)
Seamless PU upholstery Dental Chair
Auto spittoon flush and cup filler control system
Water suction and saliva ejector. (adjustable)
Import electromagnetism valve from Italy.
Control with nine memories & LCD display
Two high speed and one low speed handpieces tubes

DENTAL UNIT

BL-806



Standard Equipment

- Foot control
- Doctor stool
- Water heater system
- Luxury glass spittoon
- LED Operation Lamp with sensor
- Imported air & water tube
- Taiwan motor.
- Two three-way syringe(cold/hot water)
- Seamless PU upholstery Dental Chair
- Auto spittoon flush and cup filler control system
- Water suction and saliva ejector. (adjustable)
- Import electromagnetism valve from Italy.
- Control with nine memories & LCD display
- Two high speed and one low speed handpieces tubes

DENTAL STOOL



BL-811



BL-812



BL-815



BL-816



BL-814 A



BL-814



BL-813

DENTAL AIR COMPRESSOR



BL-901



BL-905



BL-902



BL-906

Model	Voltage (V)	Power (W)	Pressure (Mpa)	Air Flow (L/min)	Noise (dB)	Air Tank (L)	G.Weight (Kg)	Dimension (CM)	Silent Cabinet	Dryer
BL-901	AC110 /220	550	0.8	130	50	40	28	44x41x70	Options	
BL-902	AC110 /220	1100	0.8	260	52	60	53	72x40x72		
BL-903	AC110 /220	1650	0.8	390	56	90	75	97x70x40		
BL-905	AC110 /220	840	0.8	150	50	40	32	44x41x70		
BL-906	AC110 /220	1680	0.8	300	52	60	69	72x40x72		
BL-907	AC110 /220	2520	0.8	450	57	90	84	97x70x40		

DENTAL SUCTION SYSTEM



BL-702 & BL-703



BL-701



BL-704

Model	Voltage (V)	Power (W)	Max Vacuum (Kpa)	Max Pressure (Kpa)	Flow (L/min)	G.Weight (Kg)	Noise (dB)	Revolution (r/min)	Dimension (CM)	Use for Dental Unit
BL-701	AC110/220	370	-11	11.5	300	29	53	2800	39X36X86	1(100%)
BL-702	AC110/220	750	-12	12	300	36	62	2800	54X38X70	1(100%) 2(60%)
BL-703	AC110/220	1100	-18	19	600	46	62	2800	54X35X70	2(100%) 4(60%)
BL-704	AC110/220	1300	-19	19	1200	78	64	2800	100X52X76	4(100%) 6(60%)

DENTAL X-RAY



BL-101 Mobile X-ray

BL-102 Wall X-ray

Specification

Power supply: AC 110V/60Hz, AC 220V±10%, 50Hz 1kVA
 Tube voltage: 60KVP
 Tube current: 8mA
 Focus size: 0.8mm
 Total filtration: 2.5mmAL
 Exposure time: 0.1-4sec
 Leak radiation: outside one meter ≤ 0.002mGy/h
 Weight: 54Kgs
 Package Size: 139x25x52cm; 77x67x10cm
 Color: Pink, Blue, Green

Characteristic

1. Use the international advanced technology, high-efficiency, clear imaging, low leak radial, it is only 1% of national standard.
2. Light touched button, controlled with microcomputer, only press one button to select the expose parameter, fast and accurately.
3. Developing in light room, image in one minute, the doctor can diagnose conveniently for maximum.
4. Pneumatic up-down seat, more convenient and more comfortable.
5. It is for the use of the digital imaging system of mouth.

PORTABLE X-RAY



Specification

Charger Input Voltage: AC100V~240V±10%
Tube Voltage:60kV(Fixed)
Tube Current: 2mA(Fixed)
Time setting range:0.01-1.6 second(0.01 sec step)

X-ray Tube

Type: Stationary anode X-ray tube
Toshiba X-ray tube focal size:0.4mm
Target angle:12.5°C
Total filtration:1.6mm Al
(Inherent filtration:0.8mmAl, Fixed added filtration:0.8mmAl)
Display method: LCD Display
Frequency:30kHz
Rated Power: 120VA
Focal Spot To Skin Distance:110mm
Charger Input Voltage: AC100V~240V±10%
Input Voltage: DC25.2V
Supply Frequency: 50/60Hz±1%
Radiation Leakage: <0.025mGy/h
N.Weight:1.8KG
G.Weight:5Kgs
Package Size:33x28x28cm
Korea Battery:DC24V

BL-103

Specification


Charger Input Voltage: AC100V~240V±10%
Tube Voltage:60kV(Fixed)
Tube Current: 2mA(Fixed)
Time setting range:0.05~1.250 Sec with Digital Sensor
0.2~1.350 Sec with Film

X-ray Tube

Type: Stationary anode X-ray tube
Toshiba X-ray tube focal size:0.7mm
Target angle:19°C
Total filtration:1.6mm Al
Display method: LCD Display
Frequency:30kHz
Rated Power: 120VA
Focal Spot To Skin Distance:100mm
Input Voltage: DC16.8V
N.Weight:1.8KG
G.Weight:5Kgs
Package Size:33x28x28cm
Korea Battery:DC16.8V



BL-103A



Specification

Sensor External Dimensions: 37x24.5x2.5(mm)
Sensing Region:30x20mm(600mm2)
Sensor Pixel region: 1580 x1050 pixels
Sensing Pixel: 3.2 million pixels
Pixel size:20x20 um
Resolution:>20lp/mm
Life:>100,000 times
Line Long(Materials-PU):3m
Grayscale:12 bits
Connection: USB Standard



BL-105 Dental Sensor Unit

PORTABLE X-RAY

Specifications

Tube voltage: 60KV
Tube current: 1.5mA
Exposure time: 0.1~2s
Frequency: 30KHz
Tube focus: 0.3mm x 0.3mm
The distance from skin to cone: 130mm
Battery: DC14.8V 6400mAh
Rated power: 60VA
Charger input voltage: AC100V~240V ± 10%
Output voltage: DC16.8V
Ambient temperature: 5~40°C
Relative humidity: <75%RH
N.Weight:2kg
G.Weight: 4.8kg
Package size: 33x31x23cm



BL-104 Portable X-ray



BL-106
Auto X-ray Film Processor



BL-107
Manual X-ray Film Processor



BL-115
Lead Apron
0.5mmPB



BL-116
Panoramic Protective
0.5mmPB



BL-114
High-collared Apron Adult
0.5mmPB



BL-113
High-Collared Apron Children
0.5mmPB



BL-110
Protective Cap
0.5mmPB



BL-111
Protective Collar
0.5mmPB



BL-112
Lead Gloves
0.5mmPB



BL-108
Dental X-ray Film(D type)



BL-108
Dental X-ray Film(E type)



BL-109A
Sensor Position



BL-109
X-ray Film Position

DENTAL AUTOCLAVE



IIB 18L/IIB 23L



IIC 18L/IIC 23L



IIB 16L

Advantage

- ※. It adopts European B standard with 3 times pre-vacuum preceding vacuum drying. Residual humidity is less than 0.2%.
- ※. Vacuum measurement reaches -0.8bar. It sterilizes and sanitizes miscellaneous packed or unpacked medical instruments. Whether solid or A-type hollow, whether with multi-apertures or inset pipes. It efficiently tides air from every hollow area of any hollow instrument, such as handpiece, making sure that water steam reaches every corner.

Program	Temperature	Pressure	Pre-vacuum	Sterilization Time	Drying Time
NUDE	134℃	2.1bar	1	5min	18min
COTTON	134℃	2.1bar	3	5min	18min
PACKAGE	134℃	2.1bar	3	10min	18min
PLASTIC	134℃	2.1bar	1	5min	9min
NUDE	121℃	1.1bar	1	20min	18min
COTTON	121℃	1.1bar	3	20min	18min
PACKAGE	121℃	1.1bar	3	20min	18min
PLASTIC	121℃	1.1bar	1	20min	18min

Model	Voltage/Power	Chamber Size	Package Size(cm)	N.Weight (Kg)	G.Weight(Kg)
IIB 16L	220V/ 2000W	Ø230X360mm	70X54x48	45	55
IIB 18L	220V/ 2200W	Ø230X350mm	73X53x48	47	56
IIB 23L	220V/ 2200W	Ø230X450mm	74X53x48	51	60
IIC 18L	220V/1500W	Ø245X360mm	71X53X50	47	55
IIC 23L	220V/ 1700W	Ø245X450mm	77x53x50	53	62

DENTAL STERILIZER



BL-N18L-1



BL-N18L-2



BL-N18L-3



BL-40L
Hot Air Sterilizer 40L



BL-UV Sterilizer



BL-504A
Handpiece Lubricating



BL-503
Sealing Machine



BL-503A
Sealing Machine



BL-501
Water Distiller



BL-501A
Water Distiller



BL-502(5L)
Ultrasonic Cleaning



BL-502C (2.5L)
Ultrasonic Cleaning



BL-502D (6L)
Ultrasonic Cleaning



Sealing Sterilization Pouch
BL-505(5cmx200m)
BL-505A(7.5cmx200m)
BL-505B (10cmx200m)
BL-505C(15cmx200m)
BL-505D(20cmx200m)
BL-505E(90mmx260mm)

LED LIGHT CURING UNIT


Wireless LED Light Curing Unit

Specification

Working voltage for adapter: AC100-240V, 50/60HZ
Wave length: 420-480nm
Light intensity: 1100 mW/cm²-1600 mW/cm²
Package: 20×16×8cm
Weight: 0.55Kg



BL-405



Application Range: It can be used for any brand of resin material.



BL-401A Built in



BL-402 Wireless



BL-403 Wireless



BL-404 Wireless



BL-401B Built in



BL-406 Wireless



BL-406A Wireless



BL-408 Wireless

ULTRASONIC SCALER



BL-418 Ultrasonic Periodontal Treatment Device

Voltage: ~100V - 240V 50/60Hz
Input Power: 48-65VA
Scaler Amplitude: 1 μm ~ 100 μm
Endo File Amplitude: <4mm
Half Offset: 0.1N ~ 2N
Frequency: 28 ± 3kHz
PC Board Fuse: F1.6AL250V
Unit Fuse: F 5AL250V



BL-411

- ◆ Touch screen ultrasonic scale-fingertip control for the ultimate in hygiene and precise operation every-time
- ◆ Detachable Handpiece autoclave
- ◆ Automatic frequency tacking
- ◆ Scaling Perio Endo Function



BL-410

Detachable Handpiece autoclave
Automatic frequency tacking
Scaling Perio Endo Function



BL-413

LED Ultrasonic Scaler



BL-412



LED Lights

Aluminum alloy
Autoclave function
Fit to EMS Ultrasonic Scaler



BL-413




LED Lights

Aluminum alloy
Autoclave function
Fit to Satelec Ultrasonic Scaler



Endo, Amden Wrench



Aluminum Handle



Aluminum Wrench



Plastic Wrench



A Tips Holder



B Tips Holder



Scaler Tip Book

Used to remove light to medium-heavy supragingival calculus and bacterial plaque. EMS® G1 NSK® GN1 SATELEC® GD1 MECTRON® GM1	Uses flat edge to remove heavy supragingival calculus. EMS® G2 NSK® GN2 SATELEC® GD2 MECTRON® GM2	Used to remove the calculus and bacterial plaque from supragingival interdental areas. EMS® G3 NSK® GN3 SATELEC® GD3 MECTRON® GM3	Used to remove all supragingival calculus and bacterial plaque. EMS® G4 NSK® GN4 SATELEC® GD4 MECTRON® GM4	Used to remove all supragingival calculus and bacterial plaque. EMS® G5 NSK® GN5 SATELEC® GD5 MECTRON® GM5
Used to remove supragingival heavy calculus. EMS® G6 NSK® GN6 SATELEC® GD6 MECTRON® GM6	Used for the removal of bridges, crowns and posts by fracturing the cement. EMS® G7 NSK® GN7 SATELEC® GD7 MECTRON® GM7	Used for the removal of bridges, crowns and posts by fracturing the cement. EMS® G8 NSK® GN8 SATELEC® GD8 MECTRON® GM8	PERIODONTICS	
Used to level off the surface of root during periodontal flap surgery. EMS® P3D NSK® PN3D SATELEC® PD3D MECTRON® PM3D	Used to remove sub gingival calculus. EMS® P4 NSK® PN4 SATELEC® PD4 MECTRON® PM4	A powerful tool for the removal of calcification, this tip can be used for widening calcified roots in the coronal third. EMS® P4D NSK® PN4D SATELEC® PD4D MECTRON® PM4D	120° holder for U files with a shank diameter of 0.57mm, this tip is usually used for anterior root canal cleaning. EMS® E1 NSK® EN1 SATELEC® ED1 MECTRON® EM1	95° holder for U files with a shank diameter of 0.57mm, this tip is usually used for molar root canal cleaning. EMS® E2 NSK® EN2 SATELEC® ED2 MECTRON® EM2
With cooling spray disabled, this tip is used for lateral thermal condensation of gutta percha. EMS® E3 NSK® EN3 SATELEC® ED3 MECTRON® EM3	Used to remove the calcification and bad filling material of pulp cavity and overhang pulp chamber. EMS® E3D NSK® EN3D SATELEC® ED3D MECTRON® EM3D	Quickly removes obstructions and broken instruments in the mid upper thirds of large and straight canals. EMS® E4 NSK® EN4 SATELEC® ED4 MECTRON® EM4	Diamond coated for retreatments or extremely hard materials. EMS® E4D NSK® EN4D SATELEC® ED4D MECTRON® EM4D	Useful for quickly eliminating calculus and temporary cements. EMS® E5 NSK® EN5 SATELEC® ED5 MECTRON® EM5
Diamond Coated, with enhanced cutting and abrasive power laterally. Specially useful for removing very hard materials on the canal wall. EMS® E5D NSK® EN5D SATELEC® ED5D MECTRON® EM5D	Burr holder used to expand the root canal. EMS® E8 NSK® EN8 SATELEC® ED8 MECTRON® EM8	Burr holder used to expand the root canal. EMS® E9 NSK® EN9 SATELEC® ED9 MECTRON® EM9	Used for gentle treatment of the root canal of anterior teeth. EMS® E10 NSK® EN10 SATELEC® ED10 MECTRON® EM10	Used for efficient apical polishing of anterior teeth. EMS® E10D NSK® EN10D SATELEC® ED10D MECTRON® EM10D
CAVITY PREPARATION				
Used for gentle treatment of the root canal of anterior teeth. EMS® E11 NSK® EN11 SATELEC® ED11 MECTRON® EM11	Used for efficient apical polishing of anterior teeth. EMS® E11D NSK® EN11D SATELEC® ED11D MECTRON® EM11D	Ball diamond tip, for preparation of occlusal surfaces and cervical ridges. EMS® SB1 NSK® SBN1 SATELEC® SBD1 MECTRON® SBM1	Semi-spherical diamond tip, for finishing cervical interproximal areas without causing lesions to adjacent teeth. EMS® SB2 NSK® SBN2 SATELEC® SBD2 MECTRON® SBM2	Semi-spherical diamond tip, for finishing cervical interproximal areas without causing lesions to adjacent teeth. EMS® SB3 NSK® SBN3 SATELEC® SBD3 MECTRON® SBM3
Semi spherical diamond tip set 45° to the left. Allows the dentist to remove carious areas of teeth without damaging the adjacent teeth. EMS® SB1L NSK® SBN1L SATELEC® SBD1L MECTRON® SBM1L	Semi spherical diamond tip set 45° to the right. Allows the dentist to remove carious areas of teeth without damaging the adjacent teeth. EMS® SB1R NSK® SBN1R SATELEC® SBD1R MECTRON® SBM1R	Left-angled, for treating furcations, interdental calculus and very narrow spaces. EMS® P2L NSK® PN2L SATELEC® PD2L MECTRON® PM2L	Right-angled, for treating furcations, interdental calculus and very narrow spaces. EMS® P2R NSK® PN2R SATELEC® PD2R MECTRON® PM2R	Left-angled, for treating furcations, interdental calculus and very narrow spaces. EMS® P2LD NSK® PN2LD SATELEC® PD2LD MECTRON® PM2LD

BONE SURGERY

US1AUS1R to cut bones in hard to reach positions with the best efficiency. MECTRON® MODEL US1	Used for cutting bone, with greater precision and efficiency. Five sharp teeth, with lines marked, 7-8.5-10mm. MECTRON® MODEL US2	Used for highly sophisticated cutting of narrow delicate bones at the thin regions of the anatomy. Flat scalpel, sharp. MECTRON® MODEL US3	Used for periodontal osteometry, sinus membrane separation. MECTRON® MODEL US4	Used for the bone remodeling and collection of bone pieces. Scalpel with a reversed trapezoid tip, sharp on three sides. MECTRON® MODEL US5
SINUS LIFTING				
Used for highly sophisticated cutting of narrow delicate bones, cutting bone with more precision and efficiency. MECTRON® MODEL US6	Used for peeling off the nasal mucosa in the upper jaw. Flat smooth at the perimeter. MECTRON® MODEL UL1	Used for peeling off the nasal mucosa in the upper jaw, protecting the organs to the greatest degree. Horn shaped, flat smooth with dull edges. MECTRON® MODEL UL2	Used for non-wound bone cutting near soft tissue or sensitive regions. Ball end tip with diamond coating, medium size. MECTRON® MODEL UL3	Used for cutting bone valve and peeling off mucosa during lifting of maxillary sinus. Tip with diamond coating, medium size. MECTRON® MODEL UL4
IMPLANTATION				
95 Deg. used for cutting ligaments with minimal trauma when pulling out firm and solid teeth. Sharp at front side. MECTRON® MODEL UL5	Used to enlarge implantation area, avoiding soft tissue trauma. MECTRON® MODEL UI1	Used to enlarge implantation area, avoiding soft tissue trauma. MECTRON® MODEL UI2	Used for final bone implant sites preparation. MECTRON® MODEL UI7	For implantation, in plus size. Maximum 1.6mm. MECTRON® MODEL UI8
Implant tip, front size 2.8mm. MECTRON® MODEL UI9	Initial pilot osteotomy (in the maxilla). MECTRON® MODEL IM1	Pilot osteotomy in anterior region (in the maxilla). MECTRON® MODEL IM2A	To enlarge or to finalize the implant site Preparation; insert with double irrigation to avoid overheating (in the maxilla). MECTRON® MODEL IM3A	Pilot osteotomy in posterior region (in the maxilla). MECTRON® MODEL IM2P
To enlarge or to finalize the implant site Preparation; insert with double irrigation to avoid overheating (in the maxilla). MECTRON® MODEL IM3P	To optimize concentricity of implant site Preparation between $\phi 2$ and $\phi 3$ mm. MECTRON® MODEL IP2-3	To finalize the implant site preparation; insert with double irrigation to avoid overheating (in the maxilla). MECTRON® MODEL IM4P	To finalize the implant site preparation; insert with double irrigation to avoid overheating (in the maxilla). MECTRON® MODEL IM4A	To optimize concentricity of implant site Preparation between $\phi 3$ and $\phi 4$ mm. MECTRON® MODEL IP3-4
ENDODONTIC				
Used for endo treatment(diamond coated). MECTRON® MODEL UE1	Used for endo treatment(diamond coated). MECTRON® MODEL UE2	Used for endo treatment(diamond coated). MECTRON® MODEL UE3	Used for endo treatment(diamond coated). MECTRON® MODEL UE4	Used for cutting ligaments with the slightest wound when pulling out firm and solid teeth. MECTRON® MODEL EX1
Used for cutting ligaments with the slightest wound when pulling out firm and solid teeth. MECTRON® MODEL EX2	Used for cutting ligaments with the slightest wound when pulling out firm and solid teeth. MECTRON® MODEL EX3	PERIODONTAL		
Root surface scaling. Cleaning the root surface. MECTRON® MODEL UP1	Periodontal surgery; remove the tip of the broken root. Endodontic surgery, remove the inflammatory tissue. MECTRON® MODEL UP2	Used to remove calculus and clean the root surface. MECTRON® MODEL UP3		
Root surface debridement and root planning (diamond coated). MECTRON® MODEL UP4	Bone formation, root planning within the operation area. MECTRON® MODEL UP5	Bone formation, root planning within the operation area. MECTRON® MODEL UP6	Root preparation in periodontal surgery. MECTRON® MODEL OP6	Bone formation, root planning within the operation area. MECTRON® MODEL UP7

HIGH SPEED HANDPIECE



BL-03A LED Fiber Optic

Standard Head
Speed: 380. 000~420. 00rpm
Air pressure: 0. 22~0. 25Mpa
Head diameter: Ø 10. 65xH12. 6mm

Torque Head
Speed: 380. 000~420. 00rpm
Air pressure: 0. 22~0. 25Mpa
Head diameter: Ø 12. 4xH13. 6mm



Connect Direct to KAVO® MULTIflex®



BL-03 LED BL-03B LED



BL-01 Key Type BL-02 Push Type



BL-02A Color High Speed Handpiece



Dental Handpiece Oil Spray



BL-10A LED 45°



BL-10 45°



BL-10 LED 45°

LOW SPEED HANDPIECE



BL-E05

DC Brushless Micromotor
Power: 120W
Speed: 2,000-40,000rpm
Torque: Max 3.4Ncm
LED Light
Brightness: 32,000lux

Contra angle option: 16:1, 1:1, 1:1, 1:5



BL-E05M
Built in Electric Motor

Dental Implant Contra Angle



BL-06 20:1



BL-06A 20:1



BL-06B 1:5



BL-06C 4:1



BL-06D 16:1



BL-05A
Low Speed Handpiece



BL-05B
Low Speed Handpiece with
Light(Inner Water Spray) Set

Intra Oral Camera



BL-310A Holder



BL-306 17 inch Monitor

Teeth Whitening Machine



BL-1W Specifications

- Supply Voltage: AC 100-240V
- Power: 40W
- Wavelength: 460-490nm
- Density: 200-300mW/cm2
- LED tube: 8pcs
- Gross Weight: 15KG
- Packing Size: 75x75x20cm



BL-1WA Specifications

- Supply Voltage: AC 100-240 V 50-60 Hz 1.2A
- Nice new design;
- With 5 inch touch screen.
- Three LED light sources: Blue, Red, Blue&purple&Red, Blue & Purple.
- Three choices of illumination:
 1. 100% full power output
 2. 75% power output
 3. 50% power output
- Light source: 4-5 W blue LEDs & 2-3W red LEDs, 1W Purple LEDs.
- 430nm-490nm (Blue) & 620-640nm (RED) & 380nm-410nm (purple)



BL-2W Table Type Teeth Whitening Machine



BL-3W Chair Type Teeth Whitening Machine



BL-4W Whitening Kits



BL-5W Whitening Kits

DENTAL EQUIPMENT



BL-610 Dental Mobile Cart



BL-611 Dental Tray Cart



BL-607 Air Turbine



BL-606 Air Turbine



BL-90 Micro-motor (35000r/min)



BL-90 A Micro-motor (35000r/min)



BL-90 B Micro-motor (35000r/min)



BL-90 C Micro-motor (35000r/min)



BL-70 Bluishless Micromotor



BL-201A Amalgamator



BL-201 Amalgamator



BL-202 Amalgamator



BL-25B Ball Type Arch Trimmer



BL-28 Vacuum Former



BL-23A Powerful Dental Vibrator



BL-25A Model Trimmer






Shank information
 ISO 314 FG standard
 ISO 315 FG L
 ISO 316 FG XL
 ISO 204 right angle
 ISO 205 right angle L



Maximum recommended speed
 FG head diameter up to 018 300.000 rpm
 FG head diameter up to 023 200.000 rpm
 FG head diameter up to 027 120.000 rpm
 FG head diameter up to 023 30.000 rpm
 FG head diameter 025 and greater 20.000 rpm


FG Carbide Burs RA, Carbide Burs

FG/RA 1  <table><tr><th colspan="2">L</th><th colspan="11">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 001...</td><td>005</td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>FG</td><td></td><td>316 001...</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>204 001...</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td></td><td>005</td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>US No.</td><td></td><td></td><td>1/4</td><td>1/2</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td></tr></table>	L		ø 1/10mm											S	Order No	314 001...	005	006	008	010	012	014	016	018	021	023	FG		316 001...											FG XL		204 001...											RA			005	006	008	010	012	014	016	018	021	023	US No.			1/4	1/2	1	2	3	4	5	6	7	8	FG/RA 7  <table><tr><th colspan="2">L</th><th colspan="4">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 237 001</td><td>006</td><td>008</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td>313 237 001</td><td>006</td><td>008</td><td>010</td><td>012</td></tr><tr><td>FG shot</td><td></td><td>204 237 001</td><td>006</td><td>008</td><td>010</td><td>012</td></tr><tr><td>RA</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>329</td><td>330</td><td>331</td><td>332</td></tr></table>	L		ø 1/10mm				S	Order No	314 237 001	006	008	010	012	FG		313 237 001	006	008	010	012	FG shot		204 237 001	006	008	010	012	RA							US No.			329	330	331	332	FG/RA 2  <table><tr><th colspan="2">L</th><th colspan="11">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 010 001</td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>FG</td><td></td><td></td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>RA</td><td></td><td>204 010 001</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>33</td><td>34</td><td>35</td><td>36</td><td>37</td><td>38</td><td>39</td><td>41</td><td></td></tr></table>	L		ø 1/10mm											S	Order No	314 010 001	006	008	010	012	014	016	018	021	023	FG			006	008	010	012	014	016	018	021	023	RA		204 010 001										US No.			33	34	35	36	37	38	39	41	
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 001...	005	006	008	010	012	014	016	018	021	023																																																																																																																																																																										
FG		316 001...																																																																																																																																																																																				
FG XL		204 001...																																																																																																																																																																																				
RA			005	006	008	010	012	014	016	018	021	023																																																																																																																																																																										
US No.			1/4	1/2	1	2	3	4	5	6	7	8																																																																																																																																																																										
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 237 001	006	008	010	012																																																																																																																																																																																
FG		313 237 001	006	008	010	012																																																																																																																																																																																
FG shot		204 237 001	006	008	010	012																																																																																																																																																																																
RA																																																																																																																																																																																						
US No.			329	330	331	332																																																																																																																																																																																
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 010 001	006	008	010	012	014	016	018	021	023																																																																																																																																																																											
FG			006	008	010	012	014	016	018	021	023																																																																																																																																																																											
RA		204 010 001																																																																																																																																																																																				
US No.			33	34	35	36	37	38	39	41																																																																																																																																																																												
FG245  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 233 006</td><td>009</td></tr><tr><td>FG</td><td></td><td></td><td></td></tr><tr><td>FG shot</td><td></td><td>313 233 006</td><td>009</td></tr><tr><td>US No.</td><td></td><td></td><td>245</td></tr></table>	L		ø 1/10mm	S	Order No	314 233 006	009	FG				FG shot		313 233 006	009	US No.			245	FG/RA21  <table><tr><th colspan="2">L</th><th colspan="5">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 107 006</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td></tr><tr><td>FG</td><td></td><td>316 107 006</td><td>009</td><td>010</td><td>012</td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td>204 107 006</td><td>008</td><td>010</td><td>012</td><td>014</td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>56</td><td>57</td><td>58</td><td>59</td><td></td></tr></table>	L		ø 1/10mm					S	Order No	314 107 006	008	009	010	012	014	FG		316 107 006	009	010	012			FG XL								RA		204 107 006	008	010	012	014		US No.			56	57	58	59		FG21L  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 110 006</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 110 006</td><td>010</td><td>012</td></tr><tr><td>US No.</td><td></td><td></td><td>57L</td><td>58L</td></tr></table>	L		ø 1/10mm		S	Order No	314 110 006	010	012	FG					FG XL		316 110 006	010	012	US No.			57L	58L	FG/RA 31  <table><tr><th colspan="2">L</th><th colspan="7">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 107 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td><td>016</td></tr><tr><td>FG</td><td></td><td>316 107 007</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td>204 107 007</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>555</td><td>556</td><td>557</td><td>558</td><td>559</td><td>560</td></tr></table>	L		ø 1/10mm							S	Order No	314 107 007	008	009	010	012	014	016	FG		316 107 007							FG XL									RA		204 107 007	008	010	012	014	016		US No.			555	556	557	558	559	560																																			
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 233 006	009																																																																																																																																																																																			
FG																																																																																																																																																																																						
FG shot		313 233 006	009																																																																																																																																																																																			
US No.			245																																																																																																																																																																																			
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 107 006	008	009	010	012	014																																																																																																																																																																															
FG		316 107 006	009	010	012																																																																																																																																																																																	
FG XL																																																																																																																																																																																						
RA		204 107 006	008	010	012	014																																																																																																																																																																																
US No.			56	57	58	59																																																																																																																																																																																
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 110 006	010	012																																																																																																																																																																																		
FG																																																																																																																																																																																						
FG XL		316 110 006	010	012																																																																																																																																																																																		
US No.			57L	58L																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 107 007	008	009	010	012	014	016																																																																																																																																																																														
FG		316 107 007																																																																																																																																																																																				
FG XL																																																																																																																																																																																						
RA		204 107 007	008	010	012	014	016																																																																																																																																																																															
US No.			555	556	557	558	559	560																																																																																																																																																																														
FG31R  <table><tr><th colspan="2">L</th><th colspan="5">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 137 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 137 007</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1556</td><td>1557</td><td>1558</td><td>1559</td><td></td></tr></table>	L		ø 1/10mm					S	Order No	314 137 007	008	009	010	012	014	FG								FG XL		316 137 007						US No.			1556	1557	1558	1559		FG23  <table><tr><th colspan="2">L</th><th colspan="5">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 168 006</td><td>008</td><td>009</td><td>010</td><td>012</td><td>016</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 168 006</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>168</td><td>169</td><td>170</td><td>171</td><td>172</td></tr></table>	L		ø 1/10mm					S	Order No	314 168 006	008	009	010	012	016	FG								FG XL		316 168 006						US No.			168	169	170	171	172	FG207  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 150 001</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>957</td><td>958</td></tr></table>	L		ø 1/10mm		S	Order No	314 150 001	010	012	FG					US No.			957	958	FG23L  <table><tr><th colspan="2">L</th><th colspan="4">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 171 006</td><td>009</td><td>010</td><td>012</td><td>016</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 171 006</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>169L</td><td>170L</td><td>171L</td><td></td></tr></table>	L		ø 1/10mm				S	Order No	314 171 006	009	010	012	016	FG							FG XL		316 171 006					US No.			169L	170L	171L																																																	
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 137 007	008	009	010	012	014																																																																																																																																																																															
FG																																																																																																																																																																																						
FG XL		316 137 007																																																																																																																																																																																				
US No.			1556	1557	1558	1559																																																																																																																																																																																
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 168 006	008	009	010	012	016																																																																																																																																																																															
FG																																																																																																																																																																																						
FG XL		316 168 006																																																																																																																																																																																				
US No.			168	169	170	171	172																																																																																																																																																																															
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 150 001	010	012																																																																																																																																																																																		
FG																																																																																																																																																																																						
US No.			957	958																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 171 006	009	010	012	016																																																																																																																																																																																
FG																																																																																																																																																																																						
FG XL		316 171 006																																																																																																																																																																																				
US No.			169L	170L	171L																																																																																																																																																																																	
FG/RA21R  <table><tr><th colspan="2">L</th><th colspan="5">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 137 006</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td></tr><tr><td>FG</td><td></td><td>204 137 006</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td></td><td>010</td><td>012</td><td>014</td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1156</td><td>1157</td><td>1158</td><td>1159</td><td></td></tr></table>	L		ø 1/10mm					S	Order No	314 137 006	008	009	010	012	014	FG		204 137 006						RA			010	012	014			US No.			1156	1157	1158	1159		FG/RA23R  <table><tr><th colspan="2">L</th><th colspan="4">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 194 006</td><td>010</td><td>012</td><td>014</td><td>016</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td>204 194 006</td><td>010</td><td>012</td><td>014</td><td>016</td></tr><tr><td>US No.</td><td></td><td></td><td>1170</td><td>1171</td><td>1172</td><td></td></tr></table>	L		ø 1/10mm				S	Order No	314 194 006	010	012	014	016	FG							RA		204 194 006	010	012	014	016	US No.			1170	1171	1172		FG/RA33  <table><tr><th colspan="2">L</th><th colspan="7">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 168 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>016</td><td>021</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 168 007</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td>204 168 007</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>698</td><td>699</td><td>700</td><td>701</td><td>702</td><td>703</td></tr></table>	L		ø 1/10mm							S	Order No	314 168 007	008	009	010	012	016	021	FG									FG XL		316 168 007							RA		204 168 007							US No.			698	699	700	701	702	703	FG/RA33L  <table><tr><th colspan="2">L</th><th colspan="6">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 171 007</td><td>009</td><td>010</td><td>012</td><td>016</td><td>021</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>FG XL</td><td></td><td>316 171 007</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>699L</td><td>700L</td><td>701L</td><td>702L</td><td>703L</td></tr></table>	L		ø 1/10mm						S	Order No	314 171 007	009	010	012	016	021	FG								FG XL		316 171 007						US No.			699L	700L	701L	702L	703L												
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 137 006	008	009	010	012	014																																																																																																																																																																															
FG		204 137 006																																																																																																																																																																																				
RA			010	012	014																																																																																																																																																																																	
US No.			1156	1157	1158	1159																																																																																																																																																																																
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 194 006	010	012	014	016																																																																																																																																																																																
FG																																																																																																																																																																																						
RA		204 194 006	010	012	014	016																																																																																																																																																																																
US No.			1170	1171	1172																																																																																																																																																																																	
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 168 007	008	009	010	012	016	021																																																																																																																																																																														
FG																																																																																																																																																																																						
FG XL		316 168 007																																																																																																																																																																																				
RA		204 168 007																																																																																																																																																																																				
US No.			698	699	700	701	702	703																																																																																																																																																																														
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 171 007	009	010	012	016	021																																																																																																																																																																															
FG																																																																																																																																																																																						
FG XL		316 171 007																																																																																																																																																																																				
US No.			699L	700L	701L	702L	703L																																																																																																																																																																															
FG/RA33R  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 194 007</td><td>012</td><td>016</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>RA</td><td></td><td>204 194 007</td><td>012</td><td>016</td></tr><tr><td>US No.</td><td></td><td></td><td>1701</td><td>1702</td></tr></table>	L		ø 1/10mm		S	Order No	314 194 007	012	016	FG					RA		204 194 007	012	016	US No.			1701	1702	FG212L  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 184 072</td><td>014</td><td>016</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>7204</td><td>7205</td></tr></table>	L		ø 1/10mm		S	Order No	314 184 072	014	016	FG					US No.			7204	7205	FG21RX  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 137 019</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1931</td><td>1932</td></tr></table>	L		ø 1/10mm		S	Order No	314 137 019	010	012	FG					US No.			1931	1932	FG23RX  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 194 019</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm	S	Order No	314 194 019	012	FG				US No.																																																																																																									
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 194 007	012	016																																																																																																																																																																																		
FG																																																																																																																																																																																						
RA		204 194 007	012	016																																																																																																																																																																																		
US No.			1701	1702																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 184 072	014	016																																																																																																																																																																																		
FG																																																																																																																																																																																						
US No.			7204	7205																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 137 019	010	012																																																																																																																																																																																		
FG																																																																																																																																																																																						
US No.			1931	1932																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 194 019	012																																																																																																																																																																																			
FG																																																																																																																																																																																						
US No.																																																																																																																																																																																						
FG21RXG  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 194 019</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm	S	Order No	314 194 019	012	FG				US No.				FG21RXG  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 137 019</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm		S	Order No	314 137 019	010	012	FG					US No.					FG17  <table><tr><th colspan="2">L</th><th colspan="2">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 237 008</td><td>010</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1931</td><td>1932</td></tr></table>	L		ø 1/10mm		S	Order No	314 237 008	010	012	FG					US No.			1931	1932	FG36R  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>314 139 008</td><td>012</td></tr><tr><td>FG</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1958</td></tr></table>	L		ø 1/10mm	S	Order No	314 139 008	012	FG				US No.			1958																																																																																																															
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 194 019	012																																																																																																																																																																																			
FG																																																																																																																																																																																						
US No.																																																																																																																																																																																						
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 137 019	010	012																																																																																																																																																																																		
FG																																																																																																																																																																																						
US No.																																																																																																																																																																																						
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 237 008	010	012																																																																																																																																																																																		
FG																																																																																																																																																																																						
US No.			1931	1932																																																																																																																																																																																		
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	314 139 008	012																																																																																																																																																																																			
FG																																																																																																																																																																																						
US No.			1958																																																																																																																																																																																			
FG151  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>316 199 295</td><td>015</td></tr><tr><td>FG XL</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm	S	Order No	316 199 295	015	FG XL				US No.				FG152  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>316 210 295</td><td>015</td></tr><tr><td>FG L</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm	S	Order No	316 210 295	015	FG L				US No.				FG152G  <table><tr><th colspan="2">L</th><th colspan="1">ø 1/10mm</th></tr><tr><td>S</td><td>Order No</td><td>316 210 295</td><td>015</td></tr><tr><td>FG L</td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td></td></tr></table>	L		ø 1/10mm	S	Order No	316 210 295	015	FG L				US No.																																																																																																																																										
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	316 199 295	015																																																																																																																																																																																			
FG XL																																																																																																																																																																																						
US No.																																																																																																																																																																																						
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	316 210 295	015																																																																																																																																																																																			
FG L																																																																																																																																																																																						
US No.																																																																																																																																																																																						
L		ø 1/10mm																																																																																																																																																																																				
S	Order No	316 210 295	015																																																																																																																																																																																			
FG L																																																																																																																																																																																						
US No.																																																																																																																																																																																						

Burs for Laboratory

 HP1 =6		 HP31 =6		 HP21R =6																																																																																																																										
<table><tr><td>L</td><td colspan="12">ø 1/10mm</td></tr><tr><td>S</td><td>Order No</td><td>104 010 001</td><td>005</td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1/4</td><td>1/2</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td></tr></table>		L	ø 1/10mm												S	Order No	104 010 001	005	006	008	010	012	014	016	018	021	023	HP													US No.			1/4	1/2	1	2	3	4	5	6	7	8	<table><tr><td>L</td><td>3.6</td><td>3.8</td><td>4.1</td><td>4.1</td><td>4.5</td><td>4.5</td><td>4.8</td><td>5.0</td><td>5.2</td></tr><tr><td>S</td><td>Order No</td><td>104 107 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>021</td><td>023</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>555</td><td>556</td><td>557</td><td>558</td><td>559</td><td>560</td><td>561</td><td>562</td><td>563</td></tr></table>		L	3.6	3.8	4.1	4.1	4.5	4.5	4.8	5.0	5.2	S	Order No	104 107 007	008	009	010	012	014	016	018	021	023	HP												US No.			555	556	557	558	559	560	561	562	563	<table><tr><td>L</td><td colspan="4">ø 1/10mm</td></tr><tr><td>S</td><td>Order No</td><td>104 137 006</td><td>010</td><td>012</td><td>014</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1157</td><td>1158</td><td>1159</td></tr></table>		L	ø 1/10mm				S	Order No	104 137 006	010	012	014	HP						US No.			1157	1158	1159
L	ø 1/10mm																																																																																																																													
S	Order No	104 010 001	005	006	008	010	012	014	016	018	021	023																																																																																																																		
HP																																																																																																																														
US No.			1/4	1/2	1	2	3	4	5	6	7	8																																																																																																																		
L	3.6	3.8	4.1	4.1	4.5	4.5	4.8	5.0	5.2																																																																																																																					
S	Order No	104 107 007	008	009	010	012	014	016	018	021	023																																																																																																																			
HP																																																																																																																														
US No.			555	556	557	558	559	560	561	562	563																																																																																																																			
L	ø 1/10mm																																																																																																																													
S	Order No	104 137 006	010	012	014																																																																																																																									
HP																																																																																																																														
US No.			1157	1158	1159																																																																																																																									
 HP33 =6		 HP2 =6		 HP31R =6																																																																																																																										
<table><tr><td>L</td><td>3.6</td><td>3.8</td><td>4.1</td><td>4.1</td><td>4.5</td><td>4.8</td></tr><tr><td>S</td><td>Order No</td><td>104 168 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>016</td><td>021</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>698</td><td>699</td><td>700</td><td>701</td><td>702</td><td>703</td></tr></table>		L	3.6	3.8	4.1	4.1	4.5	4.8	S	Order No	104 168 007	008	009	010	012	016	021	HP									US No.			698	699	700	701	702	703	<table><tr><td>L</td><td>0.6</td><td>0.8</td><td>1.0</td><td>1.2</td><td>1.4</td><td>1.6</td><td>1.8</td><td>2.2</td></tr><tr><td>S</td><td>Order No</td><td>104 010 001</td><td>006</td><td>008</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>023</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>33</td><td>34</td><td>35</td><td>36</td><td>37</td><td>38</td><td>39</td><td>41</td></tr></table>		L	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.2	S	Order No	104 010 001	006	008	010	012	014	016	018	023	HP											US No.			33	34	35	36	37	38	39	41	<table><tr><td>L</td><td>3.6</td><td>3.8</td><td>4.1</td><td>4.1</td><td>4.5</td></tr><tr><td>S</td><td>Order No</td><td>104 137 007</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1556</td><td>1557</td><td>1558</td><td>1559</td><td></td></tr></table>		L	3.6	3.8	4.1	4.1	4.5	S	Order No	104 137 007	008	009	010	012	014	HP								US No.			1556	1557	1558	1559																
L	3.6	3.8	4.1	4.1	4.5	4.8																																																																																																																								
S	Order No	104 168 007	008	009	010	012	016	021																																																																																																																						
HP																																																																																																																														
US No.			698	699	700	701	702	703																																																																																																																						
L	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.2																																																																																																																						
S	Order No	104 010 001	006	008	010	012	014	016	018	023																																																																																																																				
HP																																																																																																																														
US No.			33	34	35	36	37	38	39	41																																																																																																																				
L	3.6	3.8	4.1	4.1	4.5																																																																																																																									
S	Order No	104 137 007	008	009	010	012	014																																																																																																																							
HP																																																																																																																														
US No.			1556	1557	1558	1559																																																																																																																								
 HP21 =6		 HP33L =6		 HP23R =6	 HP33R =6																																																																																																																									
<table><tr><td>L</td><td>3.6</td><td>3.8</td><td>4.1</td><td>4.1</td><td>4.5</td><td>4.5</td><td>4.8</td><td>5.2</td></tr><tr><td>S</td><td>Order No</td><td>104 107 006</td><td>008</td><td>009</td><td>010</td><td>012</td><td>014</td><td>016</td><td>018</td><td>023</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>56</td><td>57</td><td>58</td><td>59</td><td>60</td><td>61</td><td></td><td></td></tr></table>		L	3.6	3.8	4.1	4.1	4.5	4.5	4.8	5.2	S	Order No	104 107 006	008	009	010	012	014	016	018	023	HP											US No.			56	57	58	59	60	61			<table><tr><td>L</td><td>5.2</td><td>5.8</td><td>5.8</td><td>5.8</td><td>5.8</td></tr><tr><td>S</td><td>Order No</td><td>104 171 007</td><td>009</td><td>010</td><td>012</td><td>016</td><td>021</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>699L</td><td>700L</td><td>701L</td><td>702L</td><td>703L</td></tr></table>		L	5.2	5.8	5.8	5.8	5.8	S	Order No	104 171 007	009	010	012	016	021	HP								US No.			699L	700L	701L	702L	703L	<table><tr><td>L</td><td>4.1</td><td>4.1</td><td>4.5</td><td>4.5</td></tr><tr><td>S</td><td>Order No</td><td>104 194 006</td><td>010</td><td>012</td><td>014</td><td>016</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1170</td><td>1171</td><td>1172</td><td></td></tr></table>		L	4.1	4.1	4.5	4.5	S	Order No	104 194 006	010	012	014	016	HP							US No.			1170	1171	1172		<table><tr><td>L</td><td>4.1</td><td>4.5</td></tr><tr><td>S</td><td>Order No</td><td>104 194 007</td><td>012</td><td>016</td></tr><tr><td>HP</td><td></td><td></td><td></td><td></td></tr><tr><td>US No.</td><td></td><td></td><td>1701</td><td>1702</td></tr></table>	L	4.1	4.5	S	Order No	104 194 007	012	016	HP					US No.			1701	1702				
L	3.6	3.8	4.1	4.1	4.5	4.5	4.8	5.2																																																																																																																						
S	Order No	104 107 006	008	009	010	012	014	016	018	023																																																																																																																				
HP																																																																																																																														
US No.			56	57	58	59	60	61																																																																																																																						
L	5.2	5.8	5.8	5.8	5.8																																																																																																																									
S	Order No	104 171 007	009	010	012	016	021																																																																																																																							
HP																																																																																																																														
US No.			699L	700L	701L	702L	703L																																																																																																																							
L	4.1	4.1	4.5	4.5																																																																																																																										
S	Order No	104 194 006	010	012	014	016																																																																																																																								
HP																																																																																																																														
US No.			1170	1171	1172																																																																																																																									
L	4.1	4.5																																																																																																																												
S	Order No	104 194 007	012	016																																																																																																																										
HP																																																																																																																														
US No.			1701	1702																																																																																																																										


Diamond Trimmer Disc



HWX250				
D1/mm	254	254	230	230
D2/mm	25.4	32	22	19
Order No	1	2	3	4

MTX250		
D1/mm	254	254
D2/mm	25.4	32
Order No	1	2


DENTAL MATERIALS



Dental Cartridge Needles

Inch Type Metric Type Universal Type

1. Sharp Tri-bevel point for easy penetration
2. Sterilized by EO Gas
3. Package: Plastic contrainer per unit and 100pcs/box
4. Butt-end length:11mm,Max length:42mm
5. Materials: Stainless Steel cannulas,with ultrathin walls, lancet point,siliconized coating substantially reduces pain
6. Unit package:Label, production lot,production date,expiry date and size printed



Diameter Gauge	mm	Needle Color	Needle Length (mm)			
			0.4X11	0.4X16	0.4X21	0.4X25
27G	0.4	Blue	0.4X11	0.4X16	0.4X21	0.4X25
30G	0.3	Red	0.3X11	0.3X16	0.3X21	0.3X25

Disposable Dental Safety Syringe

1. Deedle size:27G,30G
2. Length of needle:11mm-42mm
3. Available cartridge size:1.8cc;2.2cc
4. Sterilized by EO gas; single use, self-aspirating syringe system, specially designed to prevent needle stick injuries.
5. With transparent barrel to allow visualization of aspiration.
6. Package: 100pcs needle individually packed with sterilization + 1pc autoclavable handle +6pcs rubber plug.



BL-52 Antifog Protective Glasses



BL-40 Disposable Micro Applicators



BL-41 Syringe Metal Tip
BL-41A Syringe Plastic Tip



BL-64 Dental Bibs
33x45cm,33x48cm



BL-49 Cotton Roll Divider(Rotate Type)



BL-48 Cotton Roll
8x38mm,10x38mm,12x38mm



BL-50 Protective Film



BL-51 Protective Film Divider



Endo Irrigation Needle

Irri-1
Close end tips,one side hole opening
Size:23G 27G 30G
Length:25mm

Irri-2
Close end tips,two side hole opening
Size:23G 27G 30G
Length:25mm

Irri-3
Close end tips,half side hole opening
Size:23G 27G 30G
Length:25mm



BL-44 Mixing Gun



BL-71 Mixing Tips



BL-46A Impression Tray



BL-46 Impression Tray



BL-45A Disposable High Suction



BL-45 Disposable Saliva Ejector



BL-76 Wood Poly-wedge



BL-77 Mouth Mirror



BL-73 Denture Box



BL-74 Disposable Bowl



BL-42A Disposable Propy Brushes



BL-42 Disposable Propy Cup



BL-72 Cheek Retractor



BL-75 Alcohol Lamp